


# Advancing Drug Discovery for Parkinson's Disease Through Development of HTS Assays Using Industrialized Human iPSC-derived Cell Models


Scott Schachte, Ali Fathi, Rebecca K. Fiene, Sarah Dickerson, Eugenia Jones, Coby B. Carlson, and Simon A. Hilcove  
FUJIFILM Cellular Dynamics, Inc., Madison, WI USA

## Abstract


Induced pluripotent stem cell (iPSC) technology has opened the possibility of taking somatic cells from virtually any human donor and converting them into virtually any cell type imaginable. With the help of the Parkinson's Progression Markers Initiative (PPMI), as part of The Michael J. Fox Foundation (MJFF), we generated iPSC lines from patients with Parkinson's disease (PD) carrying known risk-associated gene mutations (LRRK2 G2019S and GBA N370S) and clinical data supporting symptoms of PD. Additionally, we have engineered a SNCA A53T mutation from the isogenic iPSC donor. We then performed large-scale differentiations of these four (4) iPSC lines into biologically relevant midbrain dopaminergic neurons (i.e., iCell® DopaNeurons) to facilitate PD-focused assay development and drug screening. All cells showed similar marker expression and neuronal purity characteristic of dopaminergic neurons. We then seeded these cells into multiple assay formats, specifically investigating neuronal activity (multielectrode array), neurite outgrowth and degeneration, cell death, calcium signaling, bioenergetics and metabolism (Seahorse), alpha-synuclein accumulation, and GCase activity. We confirmed that known GBA-modulating compounds, such as cyclodextrin, modify the reduced GCase activity observed in iCell DopaNeurons containing PD-relevant mutations. Because the etiology of dopaminergic neuron cell death in PD is complex and involves multiple factors, this study shows that characterization and testing of different models in parallel is a worthwhile approach to improve understanding of disease progression and mechanism. Finally, the functional performance and consistency of iPSC-derived dopaminergic neurons demonstrate their potential use in drug screening and therapeutic validation.




## Characterization of iPSC-derived Dopaminergic Neurons



## GCase Enzyme Activity Assay




## PD-relevant iCell DopaNeurons are Highly Active on MEA




**Figure 2. Neuronal Network Activity on HD-MEA.** PD-relevant iCell DopaNeurons were co-cultured with iCell Astrocytes on 6-well high-density multi-electrode array (HD-MEA) plates from MaxWell Biosystems. Recordings were conducted on a MaxTwo system, and raster plots show the synchronous network activity detected by 1020 electrodes on Day 14. Approximately 400,000 iCell DopaNeurons and 100,000 iCell Astrocytes were plated on PEI/laminin-coated surface. These data show human iPSC-derived cells are highly functional and show spontaneous electrical activity on MEA.


## Delayed Network Bursting in LRRK2 G2019S iCell DopaNeurons




## Calcium Oscillation Assay with 3D Neurospheres




## Bioenergetic Analysis of PD iCell DopaNeurons



## Neurite Outgrowth Analysis of PD iCell DopaNeurons



## Disease-relevant LRRK2 Pharmacology



## Summary

These data demonstrate the utility of donor-derived LRRK2 + GBA and engineered A53T dopaminergic neurons across multiple assays utilized for HTS assays.

- The Parkinson's Disease Panel of iCell DopaNeurons display:
- Lot-to-lot consistency in purity and performance across multiple assays
  - Altered baseline MEA activity and calcium signaling
  - Reduced GCase activity and altered metabolic profiles

Cryopreserved dopaminergic neurons from disease-relevant backgrounds provide a biologically-relevant and reproducible system to expedite facilitate drug discovery and therapeutic development.

iPSC lines for the LRRK2 and GBA mutations are part of the Parkinson's Progression Markers Initiative (PPMI) iPSC cell bank.